







A complete range of connected infusion pumps

	Volumetric Pumps	Syringe Pumps
Infusion Pumps	 Agilia VP	 Agilia SP
Advanced Infusion Pumps	 Agilia VP MC Agilia VP MC WiFi	 Agilia SP MC Agilia SP MC WiFi
Infusion Pumps for Anaesthesia		 Agilia SP TIVA Agilia SP TIVA WiFi
Patient Controlled Analgesia (PCA)		 Agilia SP PCA WiFi

Agilia[®] SP MC Agilia[®] SP MC WiFi

Advanced syringe infusion pumps



Agilia, Vigilant and Centerium are registered trademarks by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information. These medical devices are regulated by health authorities. Please refer to the respective CE Declaration of Conformity for the applicable Medical Device Directive or Regulation. Pictures are not contractual.




On top of SP features and to meet the demand of a changing clinical environment, Agilila SP MC has been developed to be used intuitively and to promote more safety and comfort.



Agilia SP MC

Advanced syringe infusion pumps

Agilia SP MC offers a wide range of programming modes that suit all infusion protocols, from the simplest to the most advanced.



Adaptable to all protocols

- Agilia SP MC offers a wide range of programming modes: loading dose, micro-infusion, manual bolus, programmed bolus.
- The Agilia SP MC can receive up to 19 profiles with 200 drugs each.



Dose rate programming

- The pump can be programmed in dose rate with direct translation into flow rate in mL/h.

Optimal programming

- The infusion can be programmed from 0.1 up to 1200 mL/h.
- Multiple dose rate units available to match all clinical practices.

Robustness

- Push guard ensures arm protection against shocks.
- Glass-fiber material reinforces critical components protection.

Ensuring continuous infusion

- With its integrated handle, Agilia SP MC is extremely mobile, easy to carry and assemble. Up to three devices can be stacked.

