

Amika[®] Giving Sets

Enteral Pump Sets Overview
ENFit universal and ENFit


© Fresenius Kabi Deutschland GmbH. Reproduction and distribution - only with written authorization

7393097 (06/18/AC)



Enteral Pump Sets Overview ENFit universal and ENFit

Amika



	EasyBag			TwoLine EasyBag			EasyBag Mobile			VarioLine			VarioLine Comfort		Bag			Bag Mobile		
	Universal	ENFit	ENFit*	Universal	ENFit	ENFit*	Universal	ENFit	ENFit*	Universal	ENFit	ENFit*	Universal**	ENFit*	Universal	ENFit	ENFit*	Universal	ENFit	ENFit*

ART. NR.	7751907	7751900	7751917	7751910	7751903	7751994	7751999	7751905	7751916	7751909	7751902	7751919	7751998	7751904	7751908	7751956	7751914	7751913	7751906	7751915
----------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------	---------

UPPER END																				
ENPlus spike	•	•	•	•	•	•	•	•	•	•	•	•	•	•						
VarioLine bottle adapter										•	•	•	•	•				•	•	•
Attached Bag															•	•	•	•	•	•
LINE COMPONENTS																				
Drip chamber	•	•	•	•	•	•				•	•	•	•	•	•	•	•	•	•	•
Roller clamp				•	•	•							•	•	•	•	•	•	•	•
Medication port	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Cover on medication port			•			•			•			•	•	•			•			•
LOWER END																				
Adaptors for ENLock and LuerLock	•			•			•			•			•		•			•		
Cover on ENLock adaptor													•							
ENFit connector	•	•		•	•		•	•		•	•		•		•	•		•	•	
ENFit connector with cover			•			•			•			•		•			•			•
Tubing length (mm)	250	250	250	232	232	232	250	250	250	250	250	250	250	250	250	250	250	242	242	242
Length med. port to lower end (cm)	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40

Not made with DEHP. Abbrev.: ENF: ENFit / ENL: ENLock / LL: LuerLock / O: Oral. * Pump set with transparent cover on the lower-end ENFit connector and tethered cover on the ENFit medication port. ** VarioLine Comfort, ENFit universal, with transparent cover on lower-end ENLock adaptor and white cover on LuerLock adaptor for medication port.