



# Conox<sup>®</sup> 2D

## Smart Anaesthesia Monitoring

Data sheet

Two monitoring solutions in one device

Advanced digital processing algorithms for depth of anaesthesia monitoring

Both hypnotic effect and probability of response indices are measured from the same electrodes

## Conox<sup>®</sup> 2D

## Monitor main characteristics

#### Dimensions

215 x 150 x 65 mm

Weight 801 g

**Display size** 7" (154.08 x 85.92 mm)

### Battery

3.7 Volts Li-ion battery Capacity 3500 mAh Rechargeable

### Bluetooth®

5.0 Dual Mode, selectable, send only, up to 8 m

## Patient Cable length

3 m

## Indices performance

qCON index range

0 - 99

**qNOX index range** 0 - 99

EMG index range

SQI index range 0 - 100

BSR index range

0 - 100

#### Indices update rate Every 1 s

qCON index advisory signal

## Selectable, high and low limit, audio and visual

**qCON / qNOX indices update time** 

Current values correspond to EEG signal down to 20 s, i.e. displayed indices represent a moving average of 20 s of patient data.

## Indices trend graphing

Selectable, qCON, qNOX, EMG, BSR

## Signal representation

EEG mode, Spectrogram mode

## Lead off detection

Continuous

## Sensor impedance check

Automatic every 15 min or on request via panel button

Conox<sup>®</sup> is a registered trademark in the name of Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change or not available in your country; please contact your sales representative for the most updated information. These medical devices are regulated by health authorities. Please refer to their CE Declaration of Conformity for the applicable Medical Device Directive or Regulation. Please contact your Fresenius Kabi representative for more details. Pictures are not contractual.



Fresenius Kabi AG Else-Kröner-Str.1 61352 Bad Homburg GERMANY Phone: +49 (0) 6172 / 686-0 www.fresenius-kabi.com



Fresenius Vial S.A.S. Le Grand Chemin 38590 Brézins, FRANCE



## Technical specifications - EEG acquisition

**EEG resolution** 44.58 nV

## Input noise

2.9  $\mu$ V peak to valley

CMRR

J UD

Sampling frequency

## Input signal range

+/- 374 mV

## Resolution

24 bit ADC